

Atlantis® abutments implant compatibility chart



Manufacturer	Implant	Atlantis Abutment		Atlantis Crown		Atlantis CustomBase solution - Atlantis Abutment & Atlantis Crown		Atlantis Crown Abutment		Atlantis Conus Abutment	
				Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-denture
		Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium
BioHorizons	INTERNAL										
	Internal 3.5, 4.5, 5.7	x	x	x	x	x		x*	x	x	x*
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	x		x*	x	x	x*
	EXTERNAL										
	External implant 3.5	x**								x**	
	External implant 4.0, 5.0	x**	x					x		x**	
Biomet 3i	INTERNAL										
	Certain 4/3, 5/4	x	x	x	x	x		x*	x	x	x*
	Certain 3.25, 4.0, 5.0, 6.0	x	x	x	x	x		x*	x	x	x*
	Certain Prevail 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x	x	x		x*	x	x	x*
	Certain XP 4/5, 5/6	x	x	x	x	x		x	x	x	x*
	EXTERNAL										
	Miniplant 3.25	x	x					x	x	x	x
	Osseotite 3.75, 4.0, 5.0, 6.0	x	x					x	x	x	x
	Osseotite XP 4/5, 5/6	x	x					x	x	x	x
Camlog	INTERNAL										
	CONELOG 3.3	x		x		x		x		x	x
	CONELOG 3.8, 4.3	x		x		x		x		x	x
	CONELOG 5.0	x		x		x		x		x	x
	Screw-line Implant 3.3	x		x		x		x		x	x
	Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	x	x	x		x	x	x	x
Dentsply Sirona <small>OsseoSpeed implants are part of the Astra Tech Implant System.</small>	INTERNAL										
	Ankylos C/X - C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	x	x	x						x	
	Ankylos C/X - /X 3.5, 4.5, 5.5, 7.0 (indexed)	x	x	x	x	x				x	x
	OsseoSpeed EV 3.0 - Green	x		x		x	x	x			
	OsseoSpeed EV 3.6 - Purple	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed EV 4.2 - Yellow	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed EV 4.8 - Blue	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed EV 5.4 - Brown	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed Profile EV 4.2 Yellow	x	x	x		x		x		x	x
	OsseoSpeed Profile EV 4.8 Blue	x	x	x		x		x		x	x
	OsseoSpeed/OsseoSpeed TX 3.0 - Yellow	x		x		x		x			
	OsseoSpeed/OsseoSpeed TX 3.5, 4.0 - Aqua	x	x	x	x	x		x	x	x	x
	OsseoSpeed/OsseoSpeed TX 4.5, 5.0 - Lilac	x	x	x	x	x		x	x	x	x
	OsseoSpeed TX Profile 4.5, 5.0	x		x		x		x		x	x
	OsseoSpeed 4.0 - Blue (Previous)	x**								x**	
	Xive S 3.0	x		x		x		x			
	Xive S 3.4, 3.8, 4.5, 5.5	x	x	x	x	x		x	x	x	x
	Frialit 3.4, 3.8, 4.5, 5.5	x	x	x	x	x		x	x	x	x
Keystone Dental	INTERNAL										
	PrimaConnex SD 3.3, 3.5, Genesis 3.8	x	x	x	x	x		x	x	x	x
	PrimaConnex RD 4.0, 4.1, Genesis 4.5	x	x	x	x	x		x	x	x	x
	PrimaConnex WD 5.0, Genesis 5.5, 6.5	x	x	x	x	x		x	x	x	x

All trademarks and company names are the property of their respective owners.

* 2 day extra Turnaround time.

** 3 day extra Turnaround time.

Please turn over →

Manufacturer	Implant	Atlantis Abutment		Atlantis Crown		Atlantis CustomBase solution - Atlantis Abutment & Atlantis Crown		Atlantis Crown Abutment		Atlantis Conus Abutment	
				Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-enture
		Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium
MIS	INTERNAL										
	M4 NP 3.3	x		x		x		x		x	x
	M4 SP 3.75, 4.2	x	x	x	x	x		x	x	x	x
	M4 WP 5.0, 6.0	x	x	x	x	x		x	x	x	x
	SEVEN NP 3.3	x		x		x		x		x	x
	SEVEN SP 3.75, 4.2	x	x	x	x	x		x	x	x	x
	SEVEN WP 5.0, 6.0	x	x	x	x	x		x	x	x	x
Nobel Biocare	INTERNAL										
	NobelActive 3.0	x		x		x		x			
	NobelActive NP 3.5	x	x	x	x	x		x	x	x	x
	NobelActive RP 4.3, 5.0	x	x	x	x	x		x	x	x	x
	NobelActive WP 5.5	x	x	x	x	x		x	x	x	x
	NobelReplace NP 3.5	x	x	x	x	x		x	x	x	x
	NobelReplace RP 4.3	x	x	x	x	x		x	x	x	x
	NobelReplace WP 5.0	x	x	x	x	x		x	x	x	x
	NobelReplace 6.0	x	x	x	x	x		x	x	x	x
	NobelReplace Conical Connection NP 3.5	x	x	x	x	x		x	x	x	x
	NobelReplace Conical Connection RP 4.3, 5.0	x	x	x	x	x		x	x	x	x
	EXTERNAL										
	Brånemark System NP 3.3	x		x		x		x		x	x
	Brånemark System RP 3.75, 4.0	x	x	x		x		x	x	x	x
Brånemark System WP 5.0, 5.5	x	x	x		x		x	x	x		
Osstem	INTERNAL										
	Osstem TS Mini 3.5	x	x	x	x	x		x	x	x	x
	Osstem TS Regular 4.0, 4.5, 5.0, 6.0, 7.0	x	x	x	x	x		x	x	x	x
Straumann	INTERNAL										
	Bone Level 3.3 NC (Narrow CrossFit)	x	x	x	x	x		x	x	x	x
	Bone Level 4.1, 4.8 RC (Regular CrossFit)	x	x	x	x	x		x	x	x	x
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		x		x		x		x	x
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		x		x		x		x	x
	EXTERNAL										
	Standard Plus 3.3 NN (Narrow Neck)	x**								x**	
Zimmer Dental	INTERNAL										
	Tapered Screw-Vent 3.5, 4.5, 5.7	x	x	x	x	x		x	x	x	x
	Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	x		x		x		x		x	x
	Screw-Vent 3.3, 3.7, 4.5	x	x	x	x	x		x	x	x	x
	Screw-Vent 3.3, 3.7, 4.5 friction fit	x		x		x		x		x	x
	SwissPlus 4.8	x		x		x		x		x	x

All trademarks and company names are the property of their respective owners.

* 2 day extra Turnaround time.

** 3 day extra Turnaround time.